Queen Square MS Centre Team UCL Institute of Neurology, Faculty of Brain Sciences University College London

- **Team members (in alphabetical order):** Battiston Marco, Gandini Wheeler-Kingshott Claudia A. M., Grussu Francesco, Prados Ferran, Samson Rebecca S., Yiannakas Marios C.
- **Scanner:** Philips Achieva 3 Tesla
- Software version: 3.2.1
- **Submission:** Multi-slice fast field echo (FFE) scan, April 16th 2018

```
"ves";
Multi-transmit =
  Active channels =
                            "both";
                            "H1";
Nucleus =
Coil selection =
                            "SENSE-NV-16";
   element selection =
                            "HNACPC";
   connection =
                            "d";
                            "no";
Dual coil =
                                   "yes";
CLEAR =
  body tuned =
                            "no";
FOV
          AP(mm) =
                                   192;
       RL(mm) =
                            162;
       FH (mm) =
                            51;
Voxel size AP(mm) =
                                   0.5;
       RL(mm) =
                            0.5;
Slice thickness (mm) =
                                   3;
Recon voxel size (mm) =
                                   0.5;
Fold-over suppression =
                                   "no";
Reconstruction matrix =
                                   384;
SENSE =
                                   "yes";
  P reduction (RL) =
                            2;
  P os factor =
                            1;
Stacks =
                            1;
                            "parallel";
  type =
  slices =
                            17;
  slice gap =
                            "user defined";
      gap (mm) =
                            0:
  slice orientation =
                            "transverse";
                            "RL";
  fold-over direction =
  fat shift direction = "P";
Stack Offc. AP (P=+mm) =
                            -4.44073439;
       RL(L=+mm) =
                            4.00667763;
       FH (H=+mm) =
                            -3.98225141;
   Ang. AP(deg) =
       RL (deg) =
                            -6.99225521;
       FH (deg) =
                            0:
Minimum number of packages =
                                   1;
Slice scan order =
                            "interleaved";
Large table movement =
                                   "no";
PlanAlign =
                            "no";
REST slabs =
                            1;
```

```
shared =
                            "no";
                            "free";
  type =
  orientation =
                            "coronal";
  thickness (mm) =
                            80;
Rest Offc. AP (P=+mm) =
                            -104.569756;
       RL (L=+mm) =
                            -0.677741408;
       FH (H=+mm) =
                            -8.0724411;
   Ang. AP(deg) =
                            -0.0044703735;
       RL (deg) =
                            0.293235034;
       FH (deg) =
                            -0.873399973;
  power =
                            "1";
Shim Size AP(mm) =
                                   72.8470078;
       RL(mm) =
                            59.3023643;
       FH (mm) =
                            71;
   Offc. AP (P=+mm) =
                            -4.58936644;
       RL (L=+mm) =
                            3.43302441;
       FH (H=+mm) =
                            -3.96402216;
   Ang. AP(deg) =
                            0;
       RL(deg) =
                            -6.99225521;
       FH (deg) =
                            0;
                            "head first";
Patient position =
                            "supine";
    orientation =
Scan type =
                            "Imaging";
Scan mode =
                            "MS";
  technique =
                            "FFE";
                                  "no";
Contrast enhancement =
Acquisition mode =
                            "cartesian";
Fast Imaging mode =
                            "none";
Echoes =
                            1;
  partial echo =
                            "no";
  shifted echo =
                            "no";
TE =
                            "shortest";
Flip angle (deg) =
                            28;
                            "user defined";
TR =
  (ms) =
                            355;
Halfscan =
                            "yes";
                            "maximum";
Water-fat shift =
RF Shims =
                            "adaptive";
Shim =
                            "PB-volume";
ShimAlign =
                            "yes";
                            "no";
Fat suppression =
                            "no":
Water suppression =
MTC =
                                  "no";
Research prepulse =
                            "no";
Diffusion mode =
                            "no";
                            "high";
SAR mode =
B1 mode =
                            "default";
SAR Patient data =
                            "auto";
PNS mode =
                            "high";
                                   "maximum";
Gradient mode =
SofTone mode =
                                   "no";
Cardiac synchronization =
                            "no";
```

```
Heart rate > 250 bpm =
                                     "no";
Respiratory compensation = "no";
Navigator respiratory comp =
                                     "no";
Flow compensation =
                             "ves":
                             "default";
Temporal slice spacing =
fMRI echo stabilisation =
                             "no";
NSA =
                             8;
SMART =
                                     "ves":
Angio / Contrast enh. =
                                     "no";
                             "no";
Quantitative flow =
                             "no";
Manual start =
                                     "no";
Dynamic study =
Arterial Spin labeling =
                             "no";
Preparation phases =
                             "full";
Interactive F0 =
                             "no";
B0 field map =
                                     "no";
B1 field map =
                                     "no";
MIP/MPR =
                             "no";
                             " M", (3) " no";
Images =
Autoview image =
                             " no";
                             (4) " no";
Calculated images =
Reference tissue =
                             "Skeletal muscle";
Preset window contrast =
                             "intermediate";
Reconstruction mode =
                                     "real time";
Save raw data =
                                     "ves":
Hardcopy protocol =
                             "no";
Ringing filtering =
                             "rectangular";
Geometry correction =
                                     "default";
IF_info_seperator =
                             1634755923;
Total scan duration =
                             "09:54.3";
Rel. signal level (%) =
                                     100;
Act. TR/TE (ms) =
                             "355 / 9.8";
ACQ matrix M \times P =
                             "384 x 324";
ACQ \text{ voxel MPS (mm)} =
                                     "0.50 / 0.50 / 3.00";
REC voxel MPS (mm) =
                                     "0.50 / 0.50 / 3.00";
Scan percentage (%) =
                                     100;
Packages =
                             2;
Min. slice gap (mm) =
                                    0;
Optimal slices =
                             14;
Max. slices =
                             28;
Act. WFS (pix) / BW (Hz) = "4.535 / 95.7";
Min. WFS (pix) / Max. BW (Hz) = "0.658 / 660.0";
Min. TR/TE (ms) =
                             "222 / 9.8";
                             "< 55 %";
SAR / local torso =
Whole body / level =
                             "< 0.9 W/kg / normal";
                             "1.62 uT / 48 %":
B1 \text{ rms} =
PNS / level =
                             "100 % / 1st level";
Sound Pressure Level (dB) = 17.5852203;
```