

\\USER\\fMRI\_STUDIES\\DEWALD\\HALEH\_Medic\_Final\\medic\_pt4by3mm

TA: 5:52 PM: ISO Voxel size: 0.4×0.4×3.0 mmPAT: 2 Rel. SNR: 1.00 : me\_r

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	17
Dist. factor	0 %
Position	L2.2 A22.6 H56.3 mm
Orientation	T > C16.5 > S1.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	100 %
FoV read	160 mm
FoV phase	54.7 %
Slice thickness	3.0 mm
TR	872.0 ms
TE	20.0 ms
Averages	3
Concatenations	1
Filter	None
Coil elements	HC7;NC1,2

**Contrast - Common**

TR	872.0 ms
TE	20.0 ms
MTC	Off
Flip angle	28 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	3
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

**Resolution - Common**

FoV read	160 mm
FoV phase	54.7 %
Slice thickness	3.0 mm
Base resolution	384
Phase resolution	55 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36

**Resolution - iPAT**

Reference scan mode	Integrated
---------------------	------------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	17
Dist. factor	0 %
Position	L2.2 A22.6 H56.3 mm
Orientation	T > C16.5 > S1.4
Phase enc. dir.	A >> P
FoV read	160 mm
FoV phase	54.7 %
Slice thickness	3.0 mm
TR	872.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L2.2 A22.6 H56.3 mm
Orientation	T > C16.5 > S1.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L2.2 A22.6 H56.3
L	2.2 mm
A	22.6 mm
H	56.3 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	16.5
> S	1.4

**Geometry - Saturation**

Sat. region	1
Thickness	102 mm
Position	R0.5 P75.3 H17.7 mm
Orientation	C > T-13.2 > S0.4
Sat. region	2
Thickness	111 mm
Position	L0.3 A53.0 F13.5 mm
Orientation	C > T-14.3 > S0.3
Water suppr.	None
Special sat.	None

**System - Miscellaneous**

Positioning mode	ISO
Table position	H
Table position	56 mm
MSMA	S - C - T

**System - Miscellaneous**

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - All

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L2.2 A22.6 H56.3 mm
Orientation	T > C16.5 > S1.4
Rotation	0.00 deg
A >> P	88 mm
R >> L	160 mm
F >> H	51 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
--------------	----------

**System - Tx/Rx**

Frequency 1H	123.258297 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	872.0 ms
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Sequence - Part 1**

Introduction	On
Dimension	2D
Combined echoes	4

**Sequence - Part 1**

Flow comp.	Yes
Multi-slice mode	Interleaved
Bandwidth	260 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Normal
RF spoiling	On

**Sequence - Assistant**

Mode	Off
------	-----